



DMEM 低糖液体培养基

货号: BL308A

规格: 500ml

保存温度: 2-8°C, 避光保存, 有效期 12 个月

components	Molecular Weight	Concentration (mg/L)	mM
Amino Acids			
Glycine	75	30	0.4
L-Arginine hydrochloride	211	84	0.398
L-Cystine 2HCl	313	63	0.201
L-Histidine hydrochloride-H ₂ O	210	42	0.2
L-Isoleucine	131	105	0.801
L-Leucine	131	105	0.801
L-Lysine hydrochloride	183	146	0.797
L-Methionine	149	30	0.201
L-Phenylalanine	165	66	0.4
L-Serine	105	42	0.4
L-Threonine	119	95	0.798
L-Tryptophan	204	16	0.0784
L-Tyrosine disodium salt dihydrate	261	104	0.398
L-Valine	117	94	0.803
Vitamins			
Choline chloride	140	4	0.028
D-Calcium pantothenate	477	4	0.008
Folic Acid	441	4	0.009
Niacinamide	122	4	0.032
Pyridoxine hydrochloride	204	4	0.019
Riboflavin	376	0.4	0.001
Thiamine hydrochloride	337	4	0.011
i-Inositol	180	7.2	0.04
Inorganic Salts			

Calcium Chloride (CaCl ₂) (anhyd.)	111	200	1.801
Ferric Nitrate (Fe(NO ₃) ₃ ·9H ₂ O)	404	0.1	2.48E-04
Magnesium Sulfate (MgSO ₄) (anhyd.)	120	97.67	0.813
Potassium Chloride (KCl)	75	400	5.333
Sodium Chloride (NaCl)	58	6400	110.344
Sodium Phosphate monobasic (NaH ₂ PO ₄ -H ₂ O)	138	125	0.905
Other Components			
D-Glucose (Dextrose)	180	1000	5.555
Phenol Red	376.4	15	0.039
HEPES	0	0	0
Sodium Bicarbonate (NaHCO ₃)	84	3700	44.047
L-Glutamine	146	584	4
Sodium Pyruvate	110	110	1



Add: No.9 Yumin Street, Area B of the Airport Industrial Zone

Shunyi District, Beijing 101318 China

Toll Free: 400-600-4213

Website: www.biosharp.cn

Beijing Labgic

Technology Co., Ltd.

